

Insife HALOPV 5.0 Release Notes Final Version 2.0 Department: Product Management



Table of Contents

1.	Overview	3
2.	Upgradability	3
3.	Compatibility	4
4.	Changes since the last release	5
5.	Bugs resolved in HALOPV 5.0	. 18
6.	Features removed from HALOPV 5.0	. 24
7.	Known Issues & Workarounds	. 24
8.	Contact	. 26



1. Overview

Insife HALOPV 5.0 is a major release enhancing the capabilities of HALOPV 4.2 to become a multitenant solution. It contains several updates and changes, as detailed in section 4, to address some key identified bugs and issues, along with some limited new/ updated features.

2. Upgradability

HALOPV 5.0 can be upgraded from HALOPV 4.2. And it can also be deployed as a standalone full installation.



3. Compatibility

HALOPV 5.0 is compatible with the following systems:

- Argus 8.2.1.1 and Argus 8.2.3 (Please contact us for compatibility with earlier/later versions)
- ARISg 7.2 (Please contact us for compatibility with earlier/ later versions)
- Browser Compatibility

HALOPV requires a JavaScript enabled browser and supports the current and prior major release of **Google Chrome**, **Mozilla Firefox**, and **Microsoft Edge**. HALOPV leverages Oracle APEX framework and is compatible with browsers supported by it.

OBS: **Microsoft Internet Explorer** is no longer supported by HALOPV 4.0 onwards

The following Oracle Support documents should be considered for compatibility:

Oracle Software Web Browser Support Policy

https://docs.oracle.com/en/database/oracle/application-express/20.2/htmig/apex-installation-requirements.html#GUID-02BE4A34-B631-412C-8A82-EB92DABBACE0

At the time of release, Chrome was on version 119, and Edge was on version 119.



4. Changes since the last release

The following features have been updated or added in this release:

Area / Process / Group	Title	Description	Impact	Specification/ REQ ID
Aggregate Reporting	Flexibility to reuse aggregate reports in multiple records	A new enhancement has been introduced enabling the utilization of aggregate report configuration across multiple records through the addition of new field on the form 'Use linked configuration record'.	Minor Reusing existing feature. No impact on existing customers.	INSIFE- HALOPV-SS-39 (SS-39- PAGE- 143) HLP-4379
Aggregate Reporting	Aggregate Reporting configuration form updates	 Few of the fields are added and/or updated within the configuration form. They are listed as below: Form header renamed as "Data Query Configuration Form" "Case Classifications (Keywords)" renamed to "Keywords to Include" "Keywords to exclude" field is added "MedDRA SMQ" renamed to "MedDRA Term Collection to Include." "MedDRA Term Collection to Exclude" field is added Flexibility to select substance "Set Reporting Values" tab is added for the International Birth Date field, report header, title, and footer. This tab will be displayed based on certain defined criteria "PADER Configuration" and "Summary Tabulation" tabs will be displayed based on certain defined criteria 	Minor Renaming and Addition of few fields. No impact on existing customers.	INSIFE- HALOPV-SS-39 HLP-3113



Area / Process / Group	Title	Description	Impact	Specification/ REQ ID
Aggregate Reports	Aggregate report	The aggregate report ("ICSR Product-Event Line Listing") now includes a column	Minor	INSIFE-
	screen update	showing the C.1.8.1 - Case Number.	Addition of a field. No	HALOPV-SS-39
			impact on existing	
			customers.	HLP-6148
Agreements	Agreements Title	Updates have been made to display Agreements Title and Type on Agreements	Minor	*
		form page	Usability Enhancement.	HLP-4058
			No impact on existing	
			customers.	
Agreements	Tab 'Involved	Agreements workflow has been updated to include a new configuration to	Minor	*
	Organizational	display 'Involved Organizational Entries' tab.	Configuration	HLP-2891
	Entries' is now		Enhancement. No	
	Configurable		impact on existing	
			customers.	
Clinical Trials	Clinical Study Form	A new feature has been introduced for the purpose of automating the transfer of	Minor	INSIFE-
	Update	the selected Record Type from the dropdown during Clinical Study record	Automating a field	HALOPV-CLIN-
		creation in Clinical Trials - Clinical Studies, directly to the Study Phase field on the	selection. No impact on	STUD-TDS-420
		Clinical Study form (Periodicity and Status tab).	existing customers.	
				HLP-6120
Clinical Trials	Clinical Project	A new form named "Clinical Project Details" has been introduced for the Clinical	Minor	INSIFE-
	Details Form	Trials module.	New Functionality. No	HALOPV-CLIN-
	addition	This form comprises of the following sections and fields:	impact on existing	STUD-TDS-420
		Project Details tab:	customers.	
		Project ID: Free text field		HLP-3269
		Description: Free text field		
		• Studies tab: A display-only list of Studies associated with the Project		



Area / Process / Group	Title	Description	Impact	Specification/ REQ ID
Clinical Trials	Clinical Study	The "Study Project ID" field has been updated to a pop-up selector that displays	Minor	INSIFE-
	Project ID Field	the Project IDs of completed records from the Clinical Trials - Projects workflow.	Enhancing field	HALOPV-CLIN-
	Update	The field also allows users to manually enter the Project ID if required.	selection. No impact on	STUD-TDS-420
			existing customers.	HLP-3250
Clinical Trials	Addition of fields	New free text fields have been added to the case form when the report type is	Minor	INSIFE-
	for report type –	"Report from Study".	New Functionality. No	HALOPV-
	'Report from Study'	Patient ID	impact on existing	ICSRs-TDS-410
		Randomization Number	customers.	(PAGE-684)
				HLP-3246
Clinical Trials	Flexibility to select	A new enhancement has been made in case report type is "Report from Study"	Minor	INSIFE-
	study arms	and the clinical study selected has multiple arms, the user has the flexibility to	Enhancing field	HALOPV-CLIN-
		select an arm from a dropdown for a Clinical Study that's been established in the	selection. No impact on	STUD-TDS-420
		Clinical Trials workflow.	existing customers.	
				HLP-3233
Clinical Trials	Study Configuration	Study Product configuration has been updated to include ability to select –	Minor	INSIFE-
	for Product Type	Product type (Investigational Medicinal Product, Comparator, Placebo,	Configuration	HALOPV-CLIN-
		Non-investigational Medicinal Product)	Enhancement. No	STUD-TDS-420
		If a study product is blinded or not	impact on existing	
		• Study Arm(s), including designation of study products to the respective	customers.	HLP-3242 HLP-
		Study Arm(s)		3243 HLP-
				3244
Document	Attachment Page	New enhancements have been made to the document page with the	Minor	HLP-3252 *
Management	PDF enhancement	introduction of a new switch to incorporate PDF Tron URL. This functionality	Enhancement. No	HLP-5994
wanagement		allows the editor to view PDF attachments in an external online editor, do	impact on existing	1121 3334
		changes and save the change back to HaloPV.	customers.	
		ן נוומוואבי מווע למעב נווב נוומוואב שמנא נט המוטרע.	customers.	



Area / Process / Group	Title	Description	Impact	Specification/ REQ ID
E2B	R3 Export Mapping Enhancement	R3 export mapping has been updated to specifically target the fields D.10.7.1.r.3 – Continuing and D.7.1.r.6 - Family History to ensure correct order within the E2B report and accurate alignment of null flavors.	Minor Enhancement. No impact on existing customers.	INSIFE- HALOPV-ICSR- JSON-API-EXP- MAP-SS-58 HLP-6268
E2B	Mask Reporter information	The E2B R3 has been updated to mask the reporter information – Telephone when the privacy checkbox "Hide the reporter information" is selected.	Minor Usability Enhancement. No impact on existing customers.	INSIFE- HALOPV-ICSR- JSON-API-EXP- MAP-SS-58
Entities - Organisation	Tenant level enhancements	 Enhancements have been made to ensure that as a multi-tenant user: Only records associated with the user organizations will be visible, preventing the display of records from other organizations Only studies associated with the user organization will be visible to an individual user, preventing the display of studies from other organizations User will be able to view only products specifically listed for their organization, ensuring cross-tenant product visibility 	Minor No impact on existing customers. HALOPV maintain individual data segregation.	HLP-6148 INSIFE- HALOPV-FRS- 17 (URS 05.300) INSIFE- HALOPV-SS-39 (SS-39- PAGE- 128, PAGE- 166) HLP-4857 HLP-4855 HLP-4854
Entities Organisation	New Field 'Outbound File prefix'	A new field 'Outbound File Prefix' has been added to Organization, the values in this field will come from configured code list.	Minor Addition of few fields. No impact on existing customers.	INSIFE- HALOPV- ENTITIES-TDS- 414 HLP-4440



Area / Process / Group	Title	Description	Impact	Specification/ REQ ID
Factory Lock	Factory Lock on all Forms	Factory lock has been implemented on all HALO forms.	Minor Configuration Enhancement. No impact on existing customers.	* HLP-3180
Home Screen	Dashboard Updates	 The following section have been removed from the landing screen of the user: My top 5 open tasks Processing Compliance Charts 	Minor UI Enhancement. No impact on existing customers.	INSIFE- HALOPV-SS-39 (S S-39-PAGE- 1) HLP-5493
ICSRs (Medicine Incidents)	ICSR Case form	 ICSR Case Form has been updated with below changes – Add – the following fields have been added at 'Study Identification' form 'Sender(s) and Case info' section – Observe Study Type (i. e. Study Type where Reaction(s) / Event(s) were observed) Study Registration (Registration Type, Study Registration Number and Study Registration Country) Korean use (Other Studies Type) Hide the following fields have been removed from users view - 'Coded Substances in Tria'', ''Active Ingredient(s'' Text Update – text of following field has been updated – The field''Study has been unblinde'' updated to''Study is eligible for unblindin'' 	Minor Renaming and Addition of few fields. No impact on existing customers.	INSIFE- HALOPV-CLIN- STUD-TDS-420 HLP-3251



ICSRs (Medicine	Linking ICSR with	A new feature has been introduced for the purpose of auto-populating the below	Minor	INSIFE-
Incidents)	Clinical Trial	sections and fields in the ICSR using details from the selected study in the Clinical	Automating field values.	HALOPV-
	Workflow Details	Trials workflow when a Clinical Study (and optionally an Arm) is selected in the case form.	No impact on existing customers.	ICSRs-TDS-410 (PAGE-600, PAGE-640,
		 case form. For "Sender(s) and Case Info - Study Identification" section: "Sponsor Study Number (C.5.3)" is populated from "Clinical Trials - Protocol Number" "Study Name (C.5.2)" is populated from "Clinical Trials - Protocol Title" "Study Type where Reaction(s) / Event(s) were observed" is populated from "Clinical Trials - Observe Study Type" "Study Registration" details (including Registration Type, Study Registration Number, and Country) are populated from "Clinical Trials - Study Registration" "Korean use (Other Studies Type)" is populated from "Clinical Trials - Study Registration" "Korean use (Other Studies Type)" is populated from "Clinical Trials - Korean use" For "Drugs" section: "Medicinal Product Name as reported" (G.K.2.2) is populated from "Clinical Trials - Vroduct" Product MPID "Characterization of Drug Role" (G.K.1) is populated to "Suspect" by default "Investigational Product blinded" status is populated from "Clinical Trials - (Product) Blinded" Study Additional non-investigational medical products (if checkbox is marked) 	customers.	(PAGE-600, PAGE-640, PAGE660) HLP-3234 HLP-3603
		"Product" is populated from "Clinical Trials - Product"		



Area / Process /	Title	Description	Impact	Specification/
Group				REQ ID
ICSRs (Medicine	Batch Case Re-	An enhancement has been made to allow batch re-processing of cases.	Minor	INSIFE-
Incidents)	submission Feature		Usability Enhancement.	HALOPV-SS-39
	Added		No impact on existing	
			customers.	HLP-2317
ICSRs (Medicine	Alphanumeric	The Field Patient ID now supports alphanumeric characters.	Minor	INSIFE-
Incidents)	Patient ID Support		Usability Enhancement.	HALOPV-SS-39
			No impact on existing	
			customers.	HLP-5705
ICSRs (Medicine	Addition of Manual	New features have been introduced on the "Duplicate Search" page:	Minor	INSIFE-
Incidents)	Re-run Button for		Usability Enhancement.	HALOPV-SS-39
	Duplicate Search	 Allowing user to manually re-run the duplicate search 	No impact on existing	
		 Added a column 'Matching Keywords' to highlight the keywords from 	customers.	HLP-5522
		the case that matched the duplicate search criteria		
ICSRs (Medicine	Updated XML	XML encoding has been updated for GOM_GET_XML_R3 to encode the correct	Minor	INSIFE-
Incidents)	Encoding for EMA	sender organisation in the generated XML to EMA.	Enhancement. No	HALOPV-SS-39
			impact on existing	
			customers.	HLP-5440
ICSRs (Medicine	Duplicate search	The Workflow task along with the record Revision, have been incorporated on	Minor	INSIFE-
Incidents)	screen	the duplicate search screen.	Usability Enhancement.	HALOPV-SS-39
	enhancement		No impact on existing	
			customers.	HLP-6148
ICSRs (Medicine	Enhanced Medicine	A performance boost has been introduced to the search screen for Medicine	Minor	INSIFE-
Incidents)	Incident Intake	Incident Intake. It narrows down the search to only display records currently	Usability Enhancement.	HALOPV-SS-39
	Search	open in the workflow. Additionally, search results are refined to show exclusively	No impact on existing	
	Performance	the primary serious event and the primary suspect product.	customers.	HLP-5272



Area / Process /	Title	Description	Impact	Specification/
Group				REQ ID
ICSRs (Medicine	MedDRA Coding	Updates have been made to automatically populate below values for an Event	Minor	INSIFE-
Incidents)	Update	which is coded using a medically serious term (configurable list)	Automating few field	HALOPV-SS-39
		 Select "Always serious term' checkbox 	values. No impact on	
		 Select "Yes" on "Other Medically important Condition (E.1.2.3f) 	existing customers.	HLP-4928
ICSRs (Medicine	MedDRA Coding	Updates have been made to display MedDRA term text along with MedDRA	Minor	INSIFE-
Incidents)	Update	codes for the coded Verbatim.	Usability Enhancement.	HALOPV-SS-39
			No impact on existing	
			customers.	HLP-4785
ICSRs (Medicine	Reported Term	An update is implemented to the event auto-coding process. If an event verbatim	Minor	INSIFE-
Incidents)	Event Update	is available for any of the following fields then it is not replaced by the coded	Usability Enhancement.	HALOPV-SS-39
		MedDRA term:	No impact on existing	
		E.i.1.2 - Reaction / Event as Reported	customers.	HLP-5725
		G.K.7.r.1 - Indication as Reported by Primary source		
		D.9.2.r.2 - Reported Causes of Death		
		• D.9.4.r.2 - Autopsy-determined cause of Death		
ICSRs (Medicine	Null flavours on	An enhancement ensures that when the 'Reaction Outcome' is 'Unknown', then	Minor	INSIFE-
Incidents)	Event Fields	the following fields shall be populated with the Null flavor 'No Information':	Usability Enhancement.	HALOPV-
		• E.i.3.2a - Results in Death	No impact on existing	ICSRs-TDS-410
		• E.i.3.2b - Life Threatening	customers.	(PAGE-620)
		E.i.3.2c - Caused/Prolonged Hospitalisation		
		E.i.3.2d - Disabling/Incapacitating		HLP-4927
		• E.i.3.2e - Congenital Anomaly/Birth Defects		
		E.i.3.2f - Other Medically Important Condition		
ICSRs (Medicine	Updates to capture	Updates have been made to capture the Reason of Record 'deactivation' and	Minor	INSIFE-
Incidents)	Reason for record	"Re-activation", the captured reasons will be displayed in record Notes and	Addition of few fields.	HALOPV-SS-39
	'Deactivation'	Workflow history.	No impact on existing	
			customers.	HLP-2472



Area / Process / Group	Title	Description	Impact	Specification/ REQ ID
ICSRs (Medicine	Study selection in	While creating a new ICSR via wizard, users can now select following study fields	Minor	INSIFE-
Incidents)	ICSR	-	Configuration	HALOPV-SS-39
		Study Type	Enhancement. No	
		Study Arm	impact on existing	HLP-3361
		Study Product	customers.	
ICSRs (Medicine	EU specific New	An enhancement has been made to add following new EU specific fields on Drugs	Minor	INSIFE-
Incidents)	Fields Added	Tab of ICSR form to capture –	Addition of few fields.	HALOPV-ICSR-
-		INVENTED NAME	No impact on existing	JSON-API-EXP-
		SCIENTIFIC NAME	customers.	MAP-SS-584
		TRADEMARK NAME		
		STRENGTH NAME		HLP-4622
		FORM NAME		
		these fields will be available on all the applicable modules.		
ICSRs (Medicine	Organization	A new feature has been added on 'Sender Details' page, users can now search	Minor	INSIFE-
Incidents)	selection on Sender	the configured organization and copy its details into the page using the 'Copy	Usability Enhancement.	HALOPV-SS-39
	Tab	Details from Entity' button.	No impact on existing	
			customers.	HLP-4483
ICSRs (Medicine	A new tab to	A new tab has been added to capture the Sender Details, the tab is available in	Minor	INSIFE-
Incidents)	capture Sender	'Reporter and Sender Information' part of an ICSR form. It can be used to capture	Addition of new Tab. No	HALOPV-SS-39
	Details	sender's Organization and Contact details.	impact on existing	
			customers.	HLP-4244
Mailbox	Multiple email	An update ensures to remove the one-mailbox-per-workflow limit.	Minor	INSIFE-
	support		Usability Enhancement.	HALOPV-SS-39
			No impact on existing	
			customers.	HLP-6454



Area / Process /	Title	Description	Impact	Specification/
Group				REQ ID
Notifications	Enhancement to	An enhancement has been made to send the notifications only to users who	Minor	INSIFE-
	notifications	belong to the same organization as the record.	Enhancement. No	HALOPV-SS-39
			impact on existing	
			customers.	HLP-6052
Notifications	Custom Notification	Custom Notifications have been updated to send the notification to the user who	Minor	INSIFE-
		has created the record, previously it was configured to send the notification to	Usability Enhancement.	HALOPV-SS-39
		the user who has last viewed the record.	No impact on existing	(SS-06-20)
			customers.	
				HLP-3958
Products	MedDRA Synonyms	A new enhancement integrates official multilingual synonyms from the MedDRA	Minor	INSIFE-
	Enhance Search	Desktop Browser.	Usability Enhancement.	HALOPV-
	Capabilities		No impact on existing	PROD-MGMT-
			customers.	TDS-417
				HLP-5229



Area / Process /	Title	Description	Impact	Specification/
Group				REQ ID
Products	Drug Coding	An enhancement has been made on PI Browser/ Drug Browser page, users can	Minor	INSIFE-
		now use the following new tabs to select the source, browse and code	Usability Enhancement.	HALOPV-
		accordingly –	No impact on existing	PROD-MGMT- TDS-417
		Company Products	customers.	(6.1.10.,
		WHO Drug		6.1.11.,
				6.1.12.,
				6.1.13.,
				6.1.14.,
				6.1.15.,
				6.1.16.,
				6.1.17.,
				6.1.18.,
				6.1.19.,
				6.1.20.)
				/
				HLP-6115
	Searching non-	An enhancement has been made while using Drug Browser and Searching on a	Minor	INSIFE-
Products	consecutive words	non-consecutive word. The search results have been refined in such cases and a	Usability Enhancement.	HALOPV-
		hint has been added.	No impact on existing	PROD-MGMT-
			customers.	TDS-417
				HLP-4391
Products	Updates to use	Updates have been made to allow users to use Product Filters (ability to	Minor	INSIFE-
	'Product filter' and	include/exclude specific products) while applying 'Schedule on same ingredient'.	Usability Enhancement.	HALOPV-
	'Same Ingredient'		No impact on existing	PROD-MGMT-
	together		customers.	TDS-417
				HLP-2890



Area / Process / Group	Title	Description	Impact	Specification/ REQ ID
Requirements	Tab 'Requirement	Requirements intelligence module forms have been updated to include a new	Minor	*
Intelligence	verbatim and notes'	configuration to display 'Requirement verbatim and notes' tab.	Usability Enhancement.	HLP-3147
	is now Configurable		No impact on existing	
			customers.	
Risks	Risk form	The risk management tree view and the Risks child form have been updated.	Minor	*
	enhancement		Enhancement. No	HLP-5999
			impact on existing	
			customers.	
Scheduler	Updates to record	Few updates have been made to the task scheduler configuration with the ability	Minor	*
	workflow scheduler	to define:	Usability Enhancement.	HLP-5995
		Start date and stop date for the task	No impact on existing	
		Frequency of the task	customers.	
		Workflow to initiate		
System	Configuration	A new enhancement has been made to auto-calculate the due date for a case	Minor	INSIFE-
Configuration	enhancement	based on case seriousness.	Automating a field. No	HALOPV-SS-39
			impact on existing	
			customers.	HLP-5835
System	Introduction of	A new enhancement introduces three pivotal fields for configuration and	Minor	INSIFE-
Configuration	Entity and	distribution control:	Addition of few fields.	HALOPV-SS-39
	Environment-	EntitySpecific: Denotes client-specific configuration settings that are	No impact on existing	
	Specific Controls	automatically applied to all subsequent entities	customers.	HLP-6263
	with Block	EnvironmentSpecific: Represents settings that are intrinsically tied to		
	Distribution	their individual environments, discouraging replication across different		
	Functionality	environments		
		BlockDistribution: Engineered to control the progression of new		
		distribution functions, this feature can halt the transfer of specific		
		entries to succeeding environments		



Area / Process / Group	Title	Description	Impact	Specification/ REQ ID
System	Config Switch for	A new enhancement introduces config switch "Preserve original timestamps and	Minor	INSIFE-
Configuration	Preserving or	case identifier." By default, this feature is set to "Yes". When set to Yes, it retains	Usability Enhancement.	HALOPV-SS-39
	Updating	the values of the following tags from the XML provided.	No impact on existing	
	Timestamps and	C.1.2 - Date of Creation	customers.	HLP-6258
	ICSR Case Identifier	C.1.4 - Date report was first received from source		
		• C.1.5 - Date of receipt of the most recent information for this report		
		Conversely, if set to "No", the system will update these fields with the current		
		date and time upon receipt.		
System	New switch to ICSR	A new checkbox has been incorporated, granting users the option to bypass the	Minor	INSIFE-
Configuration	merge Follow-up	mandatory event/product update logic. When a Follow-Up is merged into a pre-	Usability Enhancement.	HALOPV-SS-39
	information	existing case, the system either refreshes (overwriting events/products) or	No impact on existing	
		retains the master case's original events/products, depending on the selection	customers.	HLP-6224
		made in the new checkbox.		
System	Code List	Updates have been made to ensure Code list configuration can only be Created/	Minor	INSIFE-
Configuration	Configuration	Updated or Deleted (Soft deleted) by Global Admin User, Other users can only	Configuration	HALOPV-SS-39
		View them.	Enhancement. No	(S S-39- PAGE-
			impact on existing	97, PAGE-25,
			customers.	PAGE-26)
				HLP-3932
System	Ability to store	The system configuration has been updated with the ability to store Entity ID.	Minor	INSIFE-
Configuration	Entity ID		Usability Enhancement.	HALOPV-SS-39
			No impact on existing	(SS-39- PAGE-
			customers.	4, PAGE-17,
				PAGE-19,
				PAGE-26,
				PAGE-123)
				SS-03-01
				HLP-3928 HLP-
				3934



Area / Process /	Title	Description	Impact	Specification/
Group				REQ ID
User Administration	Updates in	The 'WORKFLOW MANAGER' role has been updated to provided additional	Minor	*
	'Workflow	access of COMPLETE TASK button, after assigning the record to a user.	Configuration	HLP-1084
	Manager' Role		Enhancement. No	
			impact on existing	
			customers.	

* - This ticked does not require an update to System Specification/ Design Specification.

5. Bugs resolved in HALOPV 5.0

Area / Process / Group	Title	Description	Impact	Specification/ Bug ID
Acknowledgements	Negative Acknowledgments from MHRA	A fix ensures that negative acknowledgments received from MHRA are accurately reflected in the system.	None. Fix only.	* HLP-6035
Attachments	Attachments with invalid character	An update ensures that file names containing invalid characters are preserved as attachments, and users can download these files in readable format.	None. Functionality enhancement.	* HLP-5788
Data Collection Programs	Master Parent Title Link Fix	A fix ensures the 'Master Parent title' link works as expected for all the Views of 'Data Collection Programs'.	None. Fix only.	* HLP-1090
Data Collection Programs	Alignment of 'Type' Selection with Relevant 'Party Roles Group' in Agreements	A fix ensures that the available roles in the 'SDEA Party role groups' section are consistent with the overall Type of the SDEA selected on the main form, streamlining role selection under specific Data collection programs	None. Functionality enhancement.	* HLP-4509
Device Incidents Processing	Refined MIR_IMPORT_XML Device Data Import	MIR_IMPORT_XML package has been updated to ensure Device information is populated accurately from the imported MIR form.	None. Functionality enhancement.	* HLP-6252



Area / Process / Group	Title	Description	Impact	Specification/ Bug ID
Device Incidents Processing	Reporter Details Save	A fix ensures that when a user switches between tabs on 'Device Reporter/ Manufacturer' page the application does not loses the 'Manufacturer Information'.	None. Fix only.	* HLP-4544
Device Incidents Processing	MHRA Device Incidents	Issue in Device Incidents processing workflow have been fixed where application was displaying error message while clicking on Device Tab.	Record revisions that have corrupt data will need be handled from halo.	* HLP-3783
E2B	E2B Mapping Logic Update for D.9.3	A fix ensures that E2B mapping is updated to import information from the tag D.9.3 - patientautopsyyesno correctly into the case form.	None. Fix only.	INSIFE- HALOPV-ICSR- R3-TO-JSON- FULL- MAPPING-296 HLP-5554
E2B	E2B Mapping logic update for E.i.2.1b	An update to the E2B mapping logic ensures that E.i.2.1b - reactionmeddrallt, full hierarchy of the event term is imported in the system.	None. Fix only.	INSIFE- HALOPV-ICSR- R3-TO-JSON- FULL- MAPPING-296 HLP-5528
E2B	E2B Mapping logic update for C.1.1	A fix ensures that cases imported from PVS with an invalid C.1.1 - Sender's (Case) Safety Report Unique Identifier, are linked with a freshly created C.1.1. Furthermore, cases featuring an invalid C.1.8.1 - Worldwide Unique Case Identification Number will be associated with the valid C.1.1.	None. Functionality enhancement.	INSIFE- HALOPV-ICSR- R3-TO-JSON- FULL- MAPPING-296 HLP-6030



Area / Process / Group	Title	Description	Impact	Specification/ Bug ID
Home Screen	Quick searches integrate with User Role	A fix ensures when a user searches for a record using Quick Search, the system limits search results to logged in "user access roles".	None. Functionality enhancement.	INSIFE- HALOPV-SS-39 (SS-03-22)
				HLP-4423
ICSRs (Medicine Incidents)	Accepting decimal values for various fields	 The following fields will now be able to accept and retain decimal values: D.2.2a - Age at Time of Onset of Reaction / Event D.2.2.1a - Gestation Period When Reaction / Event Was Observed in the Foetus D.4 - Height (cm) D.10.2.2a - Age of Parent D.10.4 - Body Weight (kg) of Parent D.10.5 - Height (cm) of Parent E.i.6a - Duration of Reaction / Event G.k.2.3.r.3a - Strength G.k.4.r.6a - Duration of Drug Administration G.k.5a - Cumulative Dose to First Reaction G.k.6a - Gestation Period at Time of Exposure 	None. Functionality enhancement.	INSIFE- HALOPV-SS-39 HLP-5463
ICSRs (Medicine Incidents)	Concomitant Medication Coding for Blinded Case	A fix ensures user can code a CM (Concomitant Medication) of a Blinded case.	None. Functionality enhancement.	INSIFE- HALOPV-SS-39 HLP-4820
ICSRs (Medicine Incidents)	PI/SP Substance Processing	A fix ensures that in case only Preferred name (PN) of the substance system is entered in the case then it is copied as the Preferred Name into the Child name.	None. Fix only.	HALOPV-SS-39 HLP-4415



Area / Process / Group	Title	Description	Impact	Specification/ Bug ID
ICSRs (Medicine	Recoding for a	A fix ensures that a user will not be able to Code a Verbatim of a	None	INSIFE-
Incidents)	Completed Case	completed case, Re-Code button has been disabled in such cases.		HALOPV-SS-39
				HLP-4292
ICSRs (Medicine	Blinded Study Records	A fix ensures that a record involving a blinded study does not display	None. Fix only.	INSIFE-
Incidents)		product names within the record list.		HALOPV-SS-39
				HLP-6431
ICSRs (Medicine	Substance Entry Form	A fix ensures that on the left pane of the Substance Entry Form, the	None. Functionality	INSIFE-
Incidents)		Coded Substance is displayed. In scenarios where this is absent, the	enhancement.	HALOPV-SS-39
		Selected Substance Name (override) field is given precedence. If		
		neither of these fields is available, then the Substance/Selected		HLP-5451
		Substance Name Verbatim is presented.		
ICSRs (Medicine	Dates with leading	A fix ensures users can save dates having leading '0'(Zero). This issue	None. Functionality	INSIFE-
Incidents)	Zero	is only applicable for deployments running on APEX 21.1 or later.	enhancement.	HALOPV-ICSRs-
				TDS-410
				(PAGE- 600,
				620, 640, 641,
				650, 660, 669,
				669, 672, 675,
				678, 681, 683,
				684, 685, 333,
				336)
				HLP-4119
ICSRs (Medicine	PI Code Transfer to	A fix ensures that when a reported drug is coded using the PI Browser:	None. Fix only.	INSIFE-
Incidents)	MPID Override	the selected PI Code is transferred correctly into the MPID Override field		HALOPV-SS-39
				HLP-5663



Area / Process / Group	Title	Description	Impact	Specification/ Bug ID
Notifications	Nullified Record Notification Trigger	A fix ensures no notifications are triggered for a nullified record.	None. Fix only.	INSIFE- HALOPV-SS-39
	Prevented			
				HLP-5180
Notifications	Overdue emails sent	A functionality refinement has been executed; an overdue notification	None. Functionality	INSIFE-
	only on record	email will now be dispatched only when modifications are made to a	enhancement.	HALOPV-SS-39
	modification	record. The system will refrain from triggering the email notification		
		upon a user accessing or opening the record.		HLP-6113
Products	Product dictionary	On the product browser, updated the "Name Variation Type filter" to	None. Functionality	INSIFE-
	enhancement	blank and not "Specific Medicinal Product (MPID)".	enhancement.	HALOPV-PROD-
				MGMT-TDS-
				417
				HLP-4954
Products	Products Unspecific	A fix to ensure that unspecific terms in products can be created,	None. Fix only.	INSIFE-
	Terms	edited, nullified, De-nullified.		HALOPV-PROD-
				MGMT-TDS-
				417
				HLP-6060
Submission	Argus Integration -	HALOPV-Argus Integration (From Argus to HALOPV) updated for case	None. Functionality	*
	Case Mismatch	submission records.	enhancement.	HLP-3370
				HLP-5040
				HLP-5041
				HLP-5042
				HLP-5044



Area / Process / Group	Title	Description	Impact	Specification/ Bug ID
Submissions	New Privacy switch within the submission details form	A fix ensures that the newly added option within the Submission details form, labelled "Show patient/hide reporter", is consistently chosen both during auto-scheduled submissions and during manual submission creation. This allows patient-specific fields to be properly displayed while reporter-associated fields remain hidden.	None. Functionality enhancement.	* HLP-6224
Scheduler	Scheduler new instance creation	A fix ensures that system handles the scenario when user schedules the instance creation for a date less than current date.	None. Functionality enhancement.	INSIFE- HALOPV-SS-39 (S S-39- PAGE- 93) HLP-4450
Workflows	Automated Due Date Recalculation on Workflow Transition	A fix ensures that when a record is transferred to next workflow state, the Due Date for the record is recalculated.	None. Fix only.	* HLP-6377

* - This ticked does not require an update to System Specification/ Design Specification.



6. Features removed from HALOPV 5.0

Support for the Internet Explorer browser was removed from HALOPV 4.0 onwards.

7. Known Issues & Workarounds

Area / Process /	Title	Description	Workaround	Bug ID
Group				
CIOMS Template	CIOMS Template Manufacturer Address	At the CIOMS template, system is picking Manufacturer Address from system config. Due to this system is limited to only one manufacturer option.	While using CIOMS template, if user requires different manufacturer information based on ICSR record scope, user shall request a custom placeholder from HALOPV Support.	HLP-7107
Device Incidents	Device Intake	"Last Updated" date time in Device Intake Search is	Users can check the audit log of the record to see the	HLP-5832
Processing	Record Search	not consistent with the Audit History of the record.	last updated date.	
ICSRs (Medicine Incidents)	ICSR Coding	ICSR Recoding Trigger, which codes the uncoded terms to an Unspecific Term by the event-driven coding principle has some issues.	Users are advised to manually code such events.	HLP-6099
ICSRs (Medicine Incidents)	Test Result Value	Application displays an additional ".0" in the "Test Result (value / qualifier) (F.r.3.2)" field.	This has no impact on the actual value, it is just a representation difference.	HLP-6048
ICSRs (Medicine Incidents)	Patient Medical History	<pre><patientepisodename> is appended to <patientepisodecomments> when the former is a Targe at them to Cards</patientepisodecomments></patientepisodename></pre>	Users can manually locate the correct code from the term and populate it in the HALO structured field.	HLP-5947
		Term rather than a Code.		



Area / Process / Group	Title	Description	Workaround	Bug ID
List of Values	List of values for	For the list of values – Administration Routes, even	The deprecated values will be marked with a prefix -	HLP-6534
	Route of Admin	though the Deprecated checkbox is checked, the	(Deprecated) in the Administration Routes list of	
		value appears on the 'Duration of Drug	values.	
		Administration' field.		
Notifications	Custom Email	Custom Email Notification does not work for	Users are advised to ignore the use of " (Double	HLP-6102
	Notification	Notification's where Notification text contains	Quotes) while configuring emails notifications.	
		"(Double Quotes).		
Print	Template list under	The template list under Print option is not limited	Users are Advised to use Inline search and select the	HLP-6355
	Print	to the Workflow (Aggregate Report).	template carefully or use the 'Generate document	
			from template' option.	
Products	Supress Coding	When a product record is re-initiated, 'suppress	Users are advised to check the checkbox manually	HLP-6686
	Checkbox	coding checkbox' is not checked.	when re-initiating a product record.	
	Activation			
Products	Products Records	While searching for a product using "Records	Users are advised to use other filter criteria's while	HLP-5969
	Search	(Products)" search screen, application is not	searching for a product.	
		considering "Special Measures" as filter criteria.		
Record Batch	Reflection of batch	The Batch actions performed on too many records	Users are advised to click the "Go" button to see the	HLP-4576
Actions	Action	takes time to reflect in the Records grid as the	outcomes.	
		Activities are asynchronous.		
Workflows	Dynamic Due Date	Dynamic Due Date calculation in WF is considering	User deleted the disabled tasks.	HLP-6222
		Disabled task while calculating Due Date.		



8. Contact

Contact Insife Service Desk (<u>support@insife.com</u>) for any query related to installation or changes in this release.